

BRIOVENT

turbine driven ventilator

Adult | Pediatric | Neonatal



Made In India

Highlights:

- Adult/Paediatric / Neonatal (starting @ 2ml)
- Turbine Driven (No Compressor required)
- 15.6 Inch Colour Display with touch screen.
- Apart of conventional modes, it has PRVC, PRVC SIMV, APRV, BIPHASIC, MMV modes
- NIV with automatic leakage compensation.
- O2 Therapy (HFNC)
- Nebuliser volume compensated
- Waveforms
- Loops
- Breathe by breathe monitoring of 31 important parameters
- Personalised profiles.
- 4 hours internal battery backup
- Breathe by breathe display of Compliance and Resistance
- Int. Peep, P O.1, Insp. Hold, Exp. Hold Manoeuvre
- Lung Recruitment (optional)
- Graphical Trend/Tabular Trend/Event Log
- Low Flow Oxygen (LPO)
- Peak Inspiratory Flow upto 200 LPM (Regulated)
- Dynamic Setting limit as per safety standard – do not allow user to do wrong setting.



Versatile –

One ventilator for Adult , Pediatric , Neonatal (2ml).

Smart –

Access of Basic, Advance and Weaning Parameter at one Glance.

Advance Modes –

PRVC , APRV, BiPhasic , MMV

Wide Selection among NIV Modes

Intelligent –

Breath by breath display of Compliance and Resistance , Lung Dynamics , Lung Recruitment, Dynamic Setting limits.

**O2 Flush –**

Configurable Increment of Fio2 for a Limited period of time as Post Suction Tool.

Dynamic Setting Limits-

Dynamic Setting gives safe Operation by not allowing operator to Set Parameters out of Safety limits.

O2 Therapy –

Let your patient be treated on Non-invasive High Flow Oxygen Therapy through Nasal Cannula .

Recruitment Maneuver-

Lung Recruitment tool provides Sustained Inflation of Lungs , Aiming to Open-Up the Collapsed Alveolis,

Touch Screen and Rotary Knob Operation –

BrioVent High Definition 15.6" Bright Display gives Crystal Clear information of Vital Patient Parameters, Touch Screen Gives ease of Access along with Rotary Knob Operation.

360 Deg viewable Alarm Indication-

Top Mounted Alarm indicator with Priority based Color code gives instant alert to Operator, The location of the indicator is visible from any angle.

High Performance Noiseless Turbine –

Highly Responsive Turbine gives seamless flow deliveries from Neonatal to Adult patient Category



Easy
Fast
Intelligent



* Specifications are subject to change and may vary from those shown

www.lifelinebiz.com

L **LIFE LINE BIZ PVT. LTD.**
Adding value to life!