



# RESmart GII **BPAP** System

U-20T / Y-25T / Y-30T

An intelligent technology, featured in the RESmart GII Y/U series, it delivers automatic solution for patients with OSA and/or respiratory insufficiency. Aligned with same platform as RESmart GII BPAP System, Y/U series adapt pressure automatically to patient needs and provide a better synchronization.



U-20T / Y-25T / Y-30T

## Targeted Volume

Pressure support automatically adjusts according to target volume as to improve hypoventilation. Real time monitoring includes pressure, flow, minute ventilation(MV), tidal volume, oxygen saturation\*(SpO<sub>2</sub>) and pulse\*.

# 01

# 02

## Ventilation efficacy

Advanced leakage compensation promises adequate volume support. Ti Control, I/E Sense and Rise time guarantee better ventilation synchronization.



## Other key features

BMC+ iCode App

SpO<sub>2</sub> Kit (Optional)

GPRS/Wi-Fi Kit (Optional. Individually or combine with SpO<sub>2</sub> Kit)

# 03

\* SpO<sub>2</sub> Kit module required

## BMC Medical Co., Ltd.

Add: 5/F Main Building, No.19 Gucheng Street West, Shijingshan, Beijing 100043, China

Tel: +86-10-51663880

Fax: +86-10-51663880 ext.810

en.bmc-medical.com