



For Sleep Apnea Hypopnea
Syndrome(SAHS)



TWINPAP

COMBINATION OF CPAP & BIPAP

BIPAP NIV MASK **FEATURES:**



- ◆ Adjustable Mask With Height & Depth Support
- ◆ Optimal Convenient For Cpap/Bipap Therapy
- ◆ Special Non-sticky Velcro Soft Belt
- ◆ Adaptive Head Gear Locks
- ◆ Removable Silicon Cushion For Easy
Leak Proof & Comfortable
- ◆ Superior Quality & Skin Friendly



Technical Data



	CPAP	TWINPAP
	C-09	B-19
Mode		
<i>M</i>	CPAP	CPAP
<i>O</i>	Auto CPAP	Auto CPAP
<i>D</i>		S
<i>E</i>		T
<i>S</i>		S/T
		Auto BIPAP
Parameters		
Pressure Range	4-20 cmH2O	4-25cmH2O
Belex	2-4cmH2O	2-4cmH2O
Ramp	0-60 min	0-60 min
Esens (1-6 levels)		✓
Isens (1-6levels)		✓
Auto-leak Compensation	✓	✓
Auto-altitude Adjustment	✓	✓
Auto-On/Off	✓	✓
Pressure Titration	✓	✓
Noise Level (10cmH2O)	<30dBA	<35dBA
Heating Humidifier		
Adjustable Humidify	✓	✓
Replaceable Water Chamber	✓	✓
Data Download		
Data Details	✓	✓
Data Summary	✓	✓
Download to PC	✓	✓