



National Medical Technology

Fast - Efficient - Trusted



EMAIL ID

nationalmt17@gmail.com



PHONE NUMBER

+91-8460675102



OFFICE ADDRESS

Tirupati Complex, FF-108, Sardar Chowk, opp.
Ankur hospital, Krishnanagar, Ahmedabad,
Gujarat 380038



Clarity



BRIOVENT
turbine driven ventilator

RESPRO



COVIDIEN

positive results for life

About Us

Established as a Sole Proprietorship firm in the year 2017 at Ahmedabad (Gujarat, India), we “National Medical Technology” are a leading Wholesale Trader of a wide range of ICU, NICU, OT, GYNEC products, we supply most reliable and trusted products, we are doing hospital trunk projects. We are taking AMC, CMC,, REPAIRING , COMPENENT (BOARD) level repairing service of all hospital machines with good quality and low cost.

Why Us?

With our expertise in this domain and our vendor's firm support, we have secured a respectable position for ourselves in this domain. Different factors for which we have attained appreciation from clients spread all across the nation are:

- » Rich Vendor Base
- » High-Quality Product Range
- » Client-Centric Approach
- » Reasonable Prices
- » Ethical Business Policies
- » Large Distribution Network
- » Transparent Dealings
- » Excellent Transport & Logistic



Certificate

Clarity



RESPRO



BRIOVENT
turbine driven ventilator



Ventilator



NMT-101

Medivent Gold Ventilator

Your Search for truly Ventilator Ends Here.

Modes : VC, PC, PRVC, VS, SIMV (VC) + PS, SIMV (PC) + PS, SIMV (PRVC) + PS, 2PAP, + PS, APRV & Autoshift interactive

- 12.1 inch colour TFT screen
- PV and FV loops
- PT, FT and VT curves
- Calculation of lung mechanics
- Both pressure and flow triggering
- Quick start mode
- Automatic O₂ cell calibration
- NIV facility
- Built in Nebulization
- Alarm data log
- Servo controlled O₂ blender
- Breath to breath correction
- Apnea backup
- Emergency 100% oxygen
- Battery backup

NMT-102

Neocare Ventilator

Your Search for truly OT Products Ends Here.

Latest Model With Large LCD Screen

Flow	2 – 25 Litre per minute
Tidal Volume	6 to 300 ml
Modes of Operation	A/C, SIMV, Spontaneous, CPAP with Backup Ventilation and Standby
Mandatory Breaths	Pressure Controlled
Inspiratory Time	0.2 to 3.0 seconds
Micro controller based	
Continuous flow with time cycling mechanism	
Both pressure and volume limiting	
Pressure & Flow triggering	
Apnea backup ventilation	
Servo controlled oxygen blender	
Emergency 100% O ₂	
Patient Data and Alarm Log	
Simultaneous display of pressure and flow curve	
Calibrated pressure bar graph	
Automatic switch over to available gas input	
Automatic Calibration of Oxygen Sensor	





NMT-103

Yonker Patient Monitor

Your Search for truly ICU Products Ends Here.

12.1 inch high-resolution color TFT LCD Screen.

Thinner and lighter design, convenient and portable for user.

Freeze instant waveform for real time S-T segment analysis, pace-maker detection.

Real time S-T segment analysis, pace-maker detection.

3-Level audio/visual alarming.

Efficient resistance to interference of defibrillator and HF knife.

Long Lasting working capacity with built-in rechargeable lithium battery.

400 groups of NIBP list, 6000 seconds ECG waveform recall, 60 alarm events records recall, 7-day trend chart in storage.

Wire/wireless central monitoring system applicable to ICU/CCU/OR.

NMT-104

Smart CMS

Your Search for truly ICU Products Ends Here.


4/8/16 bed view with click of a mouse.

User friendly interface based on Windows

Central Monitoring capability upto 5 parameters - ECG, SpO2, NIBP, Respiration Rate and Temperature.

96 hours of graphical and tabular trends for all parameters with alarm event markings.

Bidirectional communication - Ease of taking NIBP readings and ECG lead selection remotely for any bed from central station.





NMT-105

Portable Monitor (Radiant P5)

Your Search for truly ICU Products Ends Here.

Size	W 255 mm, L 135 mm, H 175 mm
Weight	2 Kgs. approx.
Display	
Screen Type	Color TFT
Display Area	7" Wide
Resolution	800 x 480
Wave Form	Upto 8 waveforms can be displayed at any time Numeric Display includes HR/PR, SpO ₂ %, NiBP with MAP, Resp, Temp, Co ₂
Trends	96 Hrs. for all parameters.
Networking	2 way Central Nursing Station connectivity.
Alarms	Adjustable alarm limits
Applications	Adult, Pediatric and Neonatal.
Electrical	
Internal Battery	11.1V, 1.8 Ah, Li-ion Rechargeable battery
Charge Time	Maximum 8 hours to full capacity.
Operation	AC/DC operation
AC Mains	220 V AC±10%, 50 Hz, 1.5A F max.
Environmental	
Operating Temp.	0 to 45 °C
Operating Humidity	5 to 95% RH, non-condensing.
Storage Temp.	0 to 55 °C




NMT-106

Multipara Monitor (Spectra Gold)

OT / ICU / CCU Grade Patient Monitor

Your Search for truly ICU Products Ends Here.

- 15" TFT large bright Screen for better visibility
- 8 real time configurable waveform display
- Single knob to browse menus/sub menus
- Design facilitates transit monitoring
- Adult, pediatric and neonate application
- Big font for better viewing from a distance
- Both audio and visual alarms for vital signs
- Night mode for patient comfort and to enhance TFT life
- Easy to recognise Alarm priority system with Color Coded Visual indication, Variable pitch and Melodious sounds





NMT-107

Spectra SMART + Multipara Monitor

Your Search for truly ICU Products Ends Here.

- 12.1" TFT large bright Screen for better visibility
- 8 real time configurable waveform display
- Single knob to browse menus/sub menus
- Design facilitates transit monitoring
- Adult, pediatric and neonate application
- Big font for better viewing from a distance
- Both audio and visual alarms for vital signs
- Night mode for patient comfort and to enhance TFT life
- Easy to recognise Alarm priority system with Color Coded Visual indication, Variable pitch and Melodious sounds



NMT-108

Monitor (Spectra SLIM +) OT / ICU / CCU / Grade Patient Monitor

Your Search for truly ICU Products Ends Here.

- 10.4" TFT large bright Screen for better visibility
- 8 real time configurable waveform display
- Single knob to browse menus/sub menus
- Design facilitates transit monitoring
- Adult, pediatric and neonate application
- Big font for better viewing from a distance
- Both audio and visual alarms for vital signs
- Night mode for patient comfort and to enhance TFT life
- Easy to recognise Alarm priority system with Color Coded Visual indication, Variable pitch and Melodious sounds





NMT-109
Convenient Portable Designed For Accuracy
Your Search for truly ICU Products Ends Here.

ACCESSORIES
Nelcor™ Portable SpO₂ Patient Monitoring System, PM10N

RELIABLE RESULTS. ROBUST DATA.
The Nelcor™ portable SpO₂ patient monitoring system, PM10N, features an ergonomic shape and simple design that make it easy to use. Wherever you need SpO₂ measurements — even in difficult conditions — you can count on this pulse oximetry system to help meet the challenge.



NMT-110
Bedside SpO₂ Patient Monitoring System
Your Search for truly ICU Products Ends Here.

FEATURES AND SPECIFICATIONS	
PERFORMANCE	
Measurement Range	
SpO ₂ :	1% to 100%
Pulse rate:	20 to 250 beats per minute (bpm)
Pulse amplitude:	0.03% to 20%
ACCURACY	
Saturation (% SpO ₂ ± 1 SD)	
Adult:	70% to 100% ± 2 digits
Neonate:	70% to 100% ± 3 digits
Low perfusion:	70% to 100% ± 2 digits
Pulse rate:	20 to 250 bpm ± 3 digits
Low perfusion:	20 to 250 bpm ± 3 digits
ELECTRICAL	
Instrument	
Power requirements:	100 to 240 VAC, 50/60Hz, 45 VA
Fuse rating:	Fast-acting 2A 32VAC/DC, Fast-acting 500mA 32VAC/50DC
Battery	
Type:	Li-Ion
Battery capacity:	Minimum of 5 hours using new, fully charged battery with no alarms





NMT-111

YRINGE PUMP

Your Search for truly ICU Products Ends Here.

- Syringe Pump

- Volume Mode
- Time Mode
- Bolus

- Infusion Pump

- Drop Mode
- Time Mode
- Bolus
- 1.5 To 2 Hours Battery Backup

Features

- Auto Detect Syringe
- Easy Operate
- Alarm System



NMT-112

ECG

Your Search for truly ICU Products Ends Here.

- Single Channel ECG Machine

- Automatic & Manual Modes.



- Three Channel ECG Machine
(Contec Brand)

- Twelve Channel ECG Machine
- A4 Size Print



Features

- Up To 200 ECG Storage
- Manual, Auto & Arrhythmia Analysis Mode.
- High Resolution Colour TFT Display.
- DC/AC Power With Built-in Rechargeable Li-ion Battery.





NMT-113

DEFIBRILLATOR

Your Search for truly ICU Products Ends Here.

- Defibrillator With ECG (Biphasic)(Indian)

- Defibrillator With ECG & Printer (Biphasic)(Indian)

Features

- 7 Inches Color LCD Display
- **Waveform:** Biphasic
- **Charge Time:** Less Than 7 Seconds
- Offers Variable Energy Level Selection. (2,3,5,7,10,20,30,50,70,100,150,200 Joules)
- Adult & Pediatric Paddles Support.
- **Optional** - Printer , Monitor ,AED Mode



NMT-114

DEFIBRILLATOR

Your Search for truly ICU Products Ends Here.

- Defibrillator With Monitor & AED (Biphasic)(Indian)



- Defibrillator Monophasic / Biphasic Refurbished (Bpl , Zoll , Hp Code Master, Lp 10 ,Philips Heartstart XL)

Features

- **Battery:** 3 Nicad Batteries,
- **Battery Support :** 30 Minutes Of Pacing
- Defibrillator Waveform: 5 Msec Monophasic Pulse
- **Energy Select:** 0, 5, 10, 20, 50, 100, 200, 300, 360 Joules.





NMT-115

OXYGEN CONCENTRATOR

Your Search for truly ICU Products Ends Here.

Germany technology oxygen 5 lpm
Nosie less and USA make oxy slive use
Make- Oxytech

SPECIFICATION

- Flowrate: 0-5L/min, Purity: 93%±3%
- Outlet Pressure(Mpa): 0.04-0.07
- Sound level (dB): <50
- Power: AC220V±10%, 50/60 ±1Hz;
- Power: <480W
- Timer: accumulating timing; present timing
- Alarm: Power failure, High & Low pressure, High temperature
- Net Weight: 25kg
- Size (mm): 365 x 375 x 600

FEATURE

- PSA (Pressure Swing Absorption) technology
- Adjustable flowrate
- Large LCD display for : accumulating timing
- present timing: from 10min to 40 hrs.
- Digital: Pressure, Temperature
- Maintenance reminding
- Nebulizer

OPTIONAL

- SpO2
- Low Purity alarm: When oxygen purity is above 92%, it will give green light; When the purity is between 85% it will give
- Yellow light; when purity is below 73%, it will give red light; Maintenance alarm: Alarm for maintenance will be heard after 3000 hours of operation

NMT-116

BIPAP & C-PAP

Your Search for truly ICU Products Ends Here.

BMC BiPAP Range

Luna Bipap-T-20A Modes CPAP,S,
Auto S. (20cm)

Luna Bipap-T-25A Modes CPAP,S,
Auto S. (25cm)

Luna Bipap-T-25A Modes CPAP,S,
ST, T (25cm)





NMT-117

BIPAP CPAP

Your Search for truly ICU Products Ends Here.

CPAP & Bi-level Therapy Systems

DreamStation positive airway pressure (PAP) sleep therapy devices are designed to be as comfortable and easy to experience as sleep is intended to be. Connecting patients and care teams, DreamStation devices empower users to embrace their care with confidence, and enable care teams to practice efficient and effective patient management.



NMT-118

Bipap Machine

Your Search for truly ICU Products Ends Here.

System One features

- Advanced event detection
- Bi-Flex for comfort
- System One Resistance Control
- Oximetry module
- Encore and DirectView patient management systems
- SD card and wired, wireless, and broadband modem capable
- Digital Auto-Trak





NMT-119

CMS6000 PATIENT MONITOR



Introduction

The monitor has abundant monitoring functions and is used for monitoring of adult, pediatric and neonate. User may select different parameter configurations according to different requirements. It displays real-time data and waveform by 8" TFT color LCD, and up to 8-channel waveform and all monitoring parameters can be displayed simultaneously. 48mm thermal printer is optional. In addition, it can be connected to the central monitoring system via wired/wireless network so as to form a network monitoring system.

This device can monitor such parameters as ECG, RESP, SpO2, PR, NIBP and dual-channel TEMP, etc. It integrates parameter measurement module, display and recorder in one device to form a compact and portable equipment. At the same time, its built-in replaceable battery provides convenience for patient moving.

Main Features

- 8" TFT color LCD, multi-language interface(Simplified Chinese, Traditional Chinese, English, French, German, Turkish, Spanish, Portuguese, Italian, Dutch, Romanian, Russian, Kazakhstan, Polish, Czech).
- Fanless design, quiet, energy-saving and clean, which reduces the possibility of cross-infection.
- All-round monitor for adult, pediatric and neonate.
- With standard interface, OxyCRG, trend chart, big character interface and view bed, convenient to observe.
- Waveform, parameter color and location can be set optionally.
- Operation mode: keys and buttons.

NMT-120

CMS7000 PATIENT MONITOR



Introduction

This device can monitor such parameters as ECG, RESP, SpO2, PR, NIBP and dual-channel TEMP, etc. It integrates parameter measurement module, display and recorder in one device to form a compact and portable equipment. At the same time, its built-in replaceable battery provides convenience for patient moving.

Main Features

- 12.1" TFT color LCD, multi-language interface(Simplified Chinese, Traditional Chinese, English, French, German, Turkish, Spanish (Europe), Spanish (Mexico), Portuguese, Italian, Dutch, Romanian, Russian, Kazakh, Polish, Czech, Bulgarian, Ukrainian)
- Fanless design, quiet, energy-saving and clean, which reduces the possibility of cross-infection.
- All-round monitor for adult, pediatric and neonate.
- With standard interface, OxyCRG, trend chart, big character interface and view bed, convenient to observe.
- Waveform, parameter color and location can be set optionally.
- Operation mode: keys and buttons.
- Display up to 8-channel waveform on one screen.
- Display 7-lead ECG waveform on one screen, and with ECG step function.
- Adopt digital SpO2 technology, anti-motion and anti-ambient light interference, and measurement can be performed under the circumstance of weak filling.
- HRV analysis function.
- NIBP measurement mode: Manual/AUTO/Continuous, storage.

NMT-169

CMS8000 PATIENT MONITOR



Introduction

The monitor has abundant functions that can be used for clinical monitoring with adult, pediatric and neonate. Users may select different parameter configuration according to different requirements. The monitor, power supplied by 100-240V~,50/60Hz, adopts 12.1" color TFT LCD displaying real-time data and waveform. It can synchronously display eight-channel waveform and full monitoring parameters equipped with an optional 48mm thermal recorder. The monitor can be connected to the central monitoring system via wire or wireless network to form a network monitoring system.

This device can monitor such parameters as ECG, RESP, NIBP, SpO2, and dual-channel TEMP, etc. It integrates parameter measurement module, display and recorder in one device to form a compact and portable equipment. Its replaceable internal battery brings lot of convenient for patient moving.

Main Features

- 12.1" TFT color LCD, multi-language interface(Simplified Chinese, Traditional Chinese, English, French, German, Turkish, Spanish, Portuguese, Italian, Dutch, Romanian, Russian, Kazakh, Polish, Czech).
- Fanless design, quiet, energy-saving and clean, which reduces the possibility of cross-infection.
- All-round monitor for adult, pediatric and neonate.
- With standard interface, oxygen graph, trend graph, big character interface and view bed, convenient to observe.
- Finish all operations by keys and knobs.
- Maximum 8-channel waveform synchronous display.
- Display 7-lead ECG waveform on one screen, cascade ECG waveform display







NMT-121

BED HEAD PANEL

Your Search for truly ICU Products Ends Here

- Bed Head Panel	
- Specification	Examination Lamp, Suction Jar
Aluminium Extruded Horizontal Bed Head Panel Rail System Humidifier Vaccine Basket	• IV Pole Nurse Call System Gas Outlets As Per Requirement

	
- Air Purifier - UV Ozone Technology to Disinfect, Deodorize, Sterilize And Purify Air As Nature Does In The Atmosphere Use for Hospitals Rooms / Wards/ ICU Ozone -1000 Mg/Hr Air Changes Per Hour -2	- ECG Holter Contious Monitor for ECG Compact Design Battery Can Collect The 48 Hour ECG data Removable TF Card (Up To 2 GB)

NMT-122


PULSE OXYMETER


Your Search for truly ICU Products Ends Here.

- Hand Held Pulse Oxymeter
- SPO2 • Pulse Rate

PFT SPIROMETER

- PFT Machine
Pulmonary Function Test
- Features:**
High Resoluatun TFT Display
In-Built Thermal Printer
Storage For 2500 Patient Data
Facility To Download Data To PC


- Spirometer Features: Portable Carrying Diagnosis Volume And Flow Rate Communication With PC Easy Operation





NMT-123

EEG, EMG

Your Search for truly ICU Products Ends Here.

- EEG Machine 24 Channel
- EEG Machine 32 Channel(Motorized)

Options Available Portable EEG MACHINE:

- 64 Channel EEG Monitoring System
- 40 Channel EEG Monitoring System

Features:

- USB Interface
- No External Power Required.
- Complete EEG Data Transfer CD/DVD Hard Disc.
- Easy Connectivity With Laptop/ Desktop.



- EMG/NCV/EP System 2 Channel


- EMG/NCV/EP System 4 Channel


Features:

Portable EMG System
Real Time Data Storage & Review Facilities.
Operating System: Windows Xp/ Me/ 98.

NMT-124 TRUEBEAT 100

- 3.5" Colour TFT display with Touch Screen
- All functions accessible through keypad
- Equipped with "TISA"-True Interpretation Sensitivity Algorithm
- Detailed Interpretation package
- Personalized tracings with user's name (user editable)
- 3/4/6 Channel Print Format *
- Quantitative analysis of the tracings
- Simultaneous 12 lead acquisition
- Compact and portable
- Link to PC (Optional)
- Simple operation (Get report with single key press)
- Patient Storage memory upto 200 patient
- Excellent printing quality with its high resolution printer
- Copy Facility





1 Channel



2 Channel

NMT-126 Bipap Machine

Your Search for truly ICU Products Ends Here.


- Selectable display formats
- Splash-resistant alphanumeric keyboard and direct function keys
- Internal memory for data storage upto 200 patients*
- Automatic and manual mode operation
- 12 leads simultaneous recording : I, II, III, aVR, aVL, aVF, and V1-V6 pre-cordials
- Integrated high resolution 300 dpi thermal printer width 210mm

NMT-125

Single Channel ECG Machine

Your Search for truly ICU Products Ends Here.

- Easy paper loading
- Compact and light weight (smallest in its category)
- Integrated rechargeable battery gives more than 100 ECG tracings on full charge
- One button operation
- Attractive bag for transport ease.
- 2 Channel print format*
- Primitive interpretation*
- Auto Sensitivity*
- All functions accessible through four button keypad
- Bright LEDs indicate important functions and settings
- LCD display for parameters
- 12 Lead Acquisition
- Lead marking on paper along with speed, sensitivity and filters
- Auto and Manual operating modes
- Quantitative analysis of the tracings
- Excellent recording quality with its high resolution





NMT-127

PC ECG

Your Search for truly ICU Products Ends Here.

- Innovative high performance sampling device for desktop as well as for laptop.
- Menu based ECG system to store, analyse, average, display and print recorded ECG data.
- User friendly hardware and software.
- Simultaneous 12 Lead Acquisition.
- Compact size and lightweight.
- Import Clarity ECG data from other location.
- Manual ECG wave measurement tool.
- Compare ECG recording of one patient with its previous recording / other patient records.
- Customized colors of ECG waveforms, markers and grid lines.
- Various printing ECG formats.
- Personal configuration profiles option.
- Keyboard shortcuts for fast operations.
- Email or burn CD of ECG Data.
- Export ECG data to JPEG format.

NMT-128

SPIRO TECH

Your Search for truly ICU Products Ends Here.

- Compact and Portable hardware
- Graphical & tabular Pre/Post Test Comparison
- Predicted equation builder is provided to add/edit the equation and to evaluate the equation
- Facility to search the patient by patient Name, ID, Reference No. and Referring Doctor /referred by
- Facility to take up to 10 tests of a particular patient
- Facility to see the FVC Trends on the main screen with the result values
- Facility to print the data in different formats as per the user's need
- Facility to generate the PDF files or print formats and can be send through email as well





NMT-129

SPIRO TECH +

Your Search for truly ICU Products Ends Here.

- 3.5" Color TFT display with Touch Screen
- All functions accessible through keypad
- Built in thermal printer
- Storage for 200 patient data
- FVC, SVC and MVV Test
- Auto Interpretation
- Indian predicted equations
- USB Pen Drive Interface for test storage
- Facility to download data to PC
- Small Size light weight and easy to carry
- Mains & battery operated
- Automatic power saving mode

NMT-130

TRUE ST STRESS TEST SYSTEM

Your Search for truly ICU Products Ends Here.

Salient Features:

1. Effective GUI & aesthetics.
2. Smooth waveform, No pixelization.
3. Accurate Beat detection and Iso, j, post J point marking & results calculation.
4. Flexible display windows, user configurable with intuitive interface.
5. Flexible print pages with total preview.
6. Ergonomic treadmill.





NMT-130

INFUSA P+

Your Search for truly ICU Products Ends Here.

- Front loading with manual plunger opening syringe pump.
- Fully Micro Processor Controlled Based System
- Automatic Syringe size detection of 5,10,20,50/60 ml.
- Portable, Durable, Light Weight, use friendly and easy to operate
- Large continuous display of four lines indicates Rate Setting, Volume Infused to Patient drug name with syringe
- Any brand of Syringe can be calibrated and use
- Big LED Display of rate. Three step easy operation.
- User can set bolous program as per patient requirement online or programmable.
- Last program continue facility
- High volumetric accuracy over the entire delivery range
- Anti Bolus function active when occlusion detected. It is helps to release pressure in patient vain.
- Wide range of adjustable flow rates from 0.1 ml to 1500ml per hour

NMT-131

TELE ICU PATIENT MONITOR

Your Search for truly ICU Products Ends Here.

- Ambulatory Patient care like never befor
- All vital patient monitoring parameters at doctors desk.
- Recobro gives you Ultimate Freedom to move
- And so much more





NMT-131

Neonatal Radiant Warmer

Your Search for truly NICU Products Ends Here.

(Neonatal Radiant Warmer with detachable baby trolley)

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- Double Digital Display for ACTUAL & SET temperature
- °C/°F settings
- Heater power indicator
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL Modes
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,

NMT-132

Neonatal Intensive Care Incubator

Your Search for truly NICU Products Ends Here.

- The controller is based on Micro computerized system.
- Large bright displays are easy to monitor and for viewing from a distance.
- Digital display system provided to measure the Infant temperature. Air temperature and to set the required temperature.
- Feather touch keys provided for easy operation.
- Bar graph LED system provided to indicate the heater output.
- Fault Indications with alarm for,
 - ✦ Baby probe failure
 - ✦ Air probe failure
 - ✦ 1°C temperature high
 - ✦ 1°C temperature low
 - ✦ High temperature cut-off
 - ✦ Power failure
- Single Wall Canopy with six port holes
- Large Front door Load
- Head Up / Down facility
- Hygrometer provided to indicate the % of humidity
- IV Stand provided





NMT-133

Neonatal Double Surface Phototherapy

Your Search for truly NICU Products Ends Here.

Double surface LED Phototherapy (NDP - LED) is a new technology in Phototherapy that ensures a safe and effective therapy in neonatal hyperbilirubinaemia management. The output of this system is much higher than Fluorescent Blue Lamps and stand type phototherapy units. Therefore this proposed double surface LED phototherapy much better than conventional Phototherapy Units.

Note :

NSP - LED : Single Surface LED Phototherapy stand type

NMT-134

Neonatal Baby Trolley

Your Search for truly NICU Products Ends Here.

Baby Trolley -NBT 200 is very unique designed and it is used in Neonatal / Pediatric ward. This baby trolley can be used with Infant warmer and Phototherapy units as well. This unit has three side Acrylic Collapsible walls and trendelenburg / reverse trendelenburg system. Trolley is fixed with castors for easy mobility.





NMT-135

Neonatal Resuscitation Unit

Your Search for truly NICU Products Ends Here.

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- °C/°F settings
- Heater power indicator
- LCD Displays all fault indications
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL/PRE-WARM Modes
- APGAR Score
- Count Down Timer for 1-59 seconds
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,
- LED Phototherapy with LCD timer
- Electrical slow suction with pressure adjustment, valve & vacuum gauge indicates the pressure
- Manual Resuscitator with accessories
- O2 & Air flow meters provided with hoses

NMT-136

Neonatal Open Care System

Your Search for truly NICU Products Ends Here.

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- Double Digital Display for ACTUAL & SET temperature
- °C/°F settings
- Heater power indicator
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL Modes
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,





NMT-137

Neonatal Radiant Warmer

Your Search for truly NICU Products Ends Here.

(Neonatal Radiant Warmer only stand unit)

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- Double Digital Display for ACTUAL & SET temperature
- °C/°F settings
- Heater power indicator
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL Modes
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,

NMT-138

Neonatal Resuscitation Unit

Your Search for truly NICU Products Ends Here.

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- °C/°F settings
- Heater power indicator
- LCD Displays all fault indications
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL/PRE-WARM Modes
- APGAR Score
- Count Down Timer for 1-59 seconds
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,
- Electrical slow suction with pressure adjustment, valve & vacuum gauge indicates the pressure
- Resuscitator with complete accessories
- O2 & Air flow meters with hoses





NMT-139

Neonatal Open Care System

Your Search for truly NICU Products Ends Here.

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- Double Digital Display for ACTUAL temperature
- °C/°F settings
- Heater power indicator
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL Modes
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,
- LED Phototherapy attachment either side of the heater box

Note :

NOCP 400 : Open Care System with Metal Body Finishing



NMT-140

Neonatal Radiant Warmer

Your Search for truly NICU Products Ends Here.

(Neonatal Radiant Warmer with attached baby trolley)

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- Double Digital Display for ACTUAL & SET temperature
- °C/°F settings
- Heater power indicator
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL Modes
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,





NMT-141

OXYM AIR / OXYGEN BLENDERS

Your Search for truly NICU Products Ends Here.



Low Flow Air Oxygen Blender

Flow meter (optional) 0-15 LPM, 0-3 LPM
0-1 LPM



High/Low Flow Air Oxygen Blender

Flow meter (optional) 0-15 LPM, 0-3 LPM
0-1 LPM



**low-flow blender
with 1 lpm
cascade flow meter**

Air Oxygen blender for Heart Lung / Neonatal Care has been specially designed with 10 lpm and 1 lpm cascade flow meter at a resolution of 100ml. $\pm 3\%$ FiO_2 accuracy at any set point and Oxygen concentration control 21 to 100%. It features gas supply failure alarm and has a minimum bleed of 2.5 lpm facilitated with bleed ON/OFF facility.



**high-flow blender
with 0 - 60 lpm
cascade flow meter**



**Blender Wall Mount
bracket**



**Blender Pole Mount
bracket**

Air Oxygen flow for high flow therapy has been specially designed with 60 LPM flow meter (30 + 30 Lpm) at precise resolution 2.5 litres. $\pm 3\%$ FiO_2 accuracy at any given set point and Oxygen Concentration control 21 to 100%. It also features an gas supply failure alarm and has a minimum bleed of 3 lpm facilitated with bleed ON/OFF facility.





NMT-142

PHOTOTHERAPY FIBRE OPTIC

Your Search for truly NICU Products Ends Here.

BILLI FLEX

Blue LED with fibre optics bundle
Hi-intensity cold light
Long life LED
Wider area of therapy
Fast breakdown of hyper bilirubinemia
Adjustable Intensity



Fiber Optic Cold Light source
Long life LED lamp
Finding pneumothorax,
Artery, Hydrocephalus



Fiber Optic Cold Light source
Long life LED lamp
Battery Operated, Handheld
Finding pneumothorax,
Artery, Hydrocephalus

NMT-143

BUBBLE CPAP (OXYCAP)

Your Search for truly NICU Products Ends Here.

In-built Air/O₂ Blender & Flow Meter
In-built Reusable Pressure relief valve
In-built Digital Manometer
Oxygen fail alarm
PEEP high alarm
PEEP low alarm
PEEP in Digital display
0.1 Cm H₂O PEEP Resolution
Heater wire Humidifier
Bubble Generator
Different size Nasal Prongs &
Head Bonnet
High rigid aluminium structure





NMT-143

BUBBLE CPAP (OXYCAP)

Your Search for truly NICU Products Ends Here.

- In-built Air/O₂ Blender & Flow Meter
- Pressure relief valve
- Heater wire Humidifier
- Bubble Generator
- Different size Nasal Prongs &
- Head Bonnet
- Stainless Steel stand



NMT-144

BUBBLE CPAP WITH SS STAND

- In-built Air/O₂ Blender & Flow Meter
- In-built Reusable Pressure relief valve
- In-built Digital Manometer
- Oxygen fail alarm
- PEEP high alarm
- PEEP low alarm
- PEEP in Digital display
- 0.1 Cm H₂O PEEP Resolution
- Heater wire Humidifier
- Bubble Generator
- Different size Nasal Prongs &
- Head Bonnet
- High rigid aluminium structure



NMT-145

BUBBLE CPAP

- Dual Flow Meter for Air & Oxygen
- Pressure relief valve
- Heater wire Humidifier
- Bubble Generator
- Different size Nasal Prongs &
- Head Bonnet
- Stainless Steel stand





NMT-144

RESPIRATORY HUMIDIFIER

Your Search for truly NICU Products Ends Here.



**Servo Control
Respiratory Humidifier
i-ten technology**
Invasive & Non Invasive Mode
Flow detection
Heater Wire



Respiratory Humidifier
High/Medium/Low Setting
Heater Wire



Heater Wire Adapter for
Disposable Circuit 2pin



Heater Wire Adapter for
Reusable Circuit



8010 Temperature Probe



Heater Wire Adapter for
Disposable Circuit 3pin



8050 Temperature Sensor
with flow detection



MR 730 Temperature Probe



Dual Heater Wire Adapter
for Disposable Circuit



**Neonatal disposable
Auto feed Humidifier chamber**



**Reusable Adult
Humidifier Chamber**
Polysulphone Material
Autoclavable



**Reusable Neonatal
Humidifier Chamber**
Polysulphone Material
Autoclavable



**Pediatric / Adult disposable
Auto feed Humidifier chamber**





NMT-145

T-PIECE RESUSCITATOR

Your Search for truly NICU Products Ends Here.



Infant T Piece Resuscitator with Air/O₂ Blender & Humidifier
PIP & PEEP Control
Over Pressure Relief
Air/O₂ Blender
Heater wire Servo control
Humidifier



Infant T Piece Resuscitator with in-built air oxygen blender
PIP & PEEP Control
Over Pressure Relief
Air and oxygen inlet port
Oxygen Percentage 21% to 100%
In-built flowmeter 0 - 15 lpm



Infant T Piece Resuscitator
PIP & PEEP Control
Over Pressure Relief



Resuscitation Mask
Size: 00, 0 and 1

NMT-146

PATIENT INTERFACE

Your Search for truly NICU Products Ends Here.



Patient Interface Training Kit



Setup Guide

(Nasal Prongs 7 sizes, Nasal Mask 4 sizes, flow tube with pillow Head measuring tape, Head bonnet 8 sizes, Interface setup Guide.)



Bubble CPAP combo kit
Contains single heater wire circuit, Disposable auto feed humidifier chamber, preset pressure manifold and reusable bubble jar



Preset Pressure Manifold



Reusable Bubble Jar



Highflow oxygen therapy breathing tube



Highflow oxygen therapy cannula with temperature probe port
Premature, Neonatal, Infant and Peadiatric sizes are available



15mm male

Highflow oxygen therapy cannula for ventilator

Premature, Neonatal, Infant and Peadiatric sizes are available






NMT-146

LED OPERATION THEATER LIGHT (HOSPI CT 3 + 3 HEX)

Your Search for truly OT Products Ends Here.

Technical Specification

Lux Power Intensity	: 1,00,000 + 1,00,000
Color Rendering Index	: 95
Number of led	: 57 + 57
Depth Lumination	: 10 - 15mm
Colour Temp.	: 4300-5200K
Dimming	: 10% ---100%
LED Average Life	: >50000 Hrs
Power Supply	: 220V/50 Drive by 24 volt
Focusing	: Adjustable By Handle



NMT-147

LED OPERATION THEATER LIGHT (HOSPI CT 3 ROUND)

Your Search for truly OT Products Ends Here.

Technical Specification

Lux Power Intensity	: 1,00,000
Color Rendering Index	: 95
Number of led	: 57
Depth Lumination	: 10 - 15mm
Colour Temp.	: 4300-5200K
Dimming	: 10% ---100%
LED Average Life	: >50000 Hrs
Power Supply	: 220V/50 Drive by 24 volt
Focusing	: Adjustable By Handle





NMT-148

LED OPERATION THEATER LIGHT (HOSPI CT 3 APPLE)

Your Search for truly OT Products Ends Here.



Technical Specification

Lux Power Intensity	: 1,10,000
Color Rendering Index	: 91
Number of led	: 54
Depth illumination	: 10 - 15mm
Colour Temp.	: 3800-5200K
Dimming	: 10%---100%
LED Average Life	: >50000 Hrs
Power Supply	: 110V/220 Drive by 24 volt
Focusing	: Adjustable By Handle

NMT-149

LED OPERATION THEATER LIGHT (HOSPI CT 3 APPLE)

Your Search for truly OT Products Ends Here.

Technical Specification

Lux Power Intensity	: 1,10,000
Color Rendering Index	: 91
Number of led	: 54
Depth illumination	: 10 - 15mm
Colour Temp.	: 3800-5200K
Dimming	: 10%---100%
LED Average Life	: >50000 Hrs
Power Supply	: 110V/220 Drive by 24 volt
Focusing	: Adjustable By Handle





NMT-150

LED OPERATION THEATER LIGHT (HOSPI CT-5 BUTTERFLY)

Your Search for truly OT Products Ends Here.



Technical Specification

Lux Power Intensity	: 1,60,000
Color Rendering Index	: 91
Number of led	: 105
Field Diameter	: 125mm
Colour Temp.	: 3800-4500K
Dimming	: 10%---100%
LED Average Life	: >50000 Hrs
Power Supply	: 110V/220 Drive by 24 volt
Focusing	: Adjustable By Handle

LED MAKE **OSRAM**


NMT-151

LED OPERATION THEATER LIGHT (HOSPI CT 4 + 4 HEX)

Your Search for truly OT Products Ends Here.

Technical Specification

Lux Power Intensity	: 1,20,000 + 1,20,000
Color Rendering Index	: 95
Number of led	: 76 + 76
Depth Lumination	: 10 - 15mm
Colour Temp.	: 4300-5200K
Dimming	: 10% ---100%
LED Average Life	: >50000 Hrs
Power Supply	: 220V/50 Drive by 24 volt
Focusing	: Adjustable By Handle



**NMT-152****LED OPERATION THEATER LIGHT (HOSPI CT-4 BUTTERFLY)**

Your Search for truly OT Products Ends Here.

Technical Specification


Lux Power Intensity	: 1,30,000
Color Rendering Index	: 91
Number of led	: 84
Field Diameter	: 125mm
Colour Temp.	: 3800-4500K
Dimming	: 10%---100%
LED Average Life	: >50000 Hrs
Power Supply	: 110V/220 Drive by 24 volt
Focusing	: Adjustable By Handle

**NMT-153****LED OPERATION THEATER LIGHT (HOSPI CT M5 HEX)**

Your Search for truly OT Products Ends Here.

Technical Specification

Lux Power Intensity	: 1,60,000 + 1,60,000
Color Rendering Index	: 95
Number of led	: 95 + 95
Depth Lumination	: 10 - 15mm
Colour Temp.	: 4300-5200K
Dimming	: 10% ---100%
LED Average Life	: >50000 Hrs
Power Supply	: 220V/50 Drive by 24 volt
Focusing	: Remote Operated





NMT-154

LED OPERATION THEATER LIGHT (HOSPI CT-4 BUTTERFLY)

Your Search for truly OT Products Ends Here.



Technical Specification

Lux Power Intensity	: 1,30,000
Color Rendering Index	: 91
Number of led	: 84
Field Diameter	: 125mm
Colour Temp.	: 3800-4500K
Dimming	: 10%---100%
LED Average Life	: >50000 Hrs
Power Supply	: 110V/220 Drive by 24 volt
Focusing	: Adjustable By Handle

NMT-155

MEDIVENT EXCEL

Your Search for truly OT Products Ends Here.

- Micro-controller based
- No Need of driving gas
- Apnea backup ventilation
- Adjustable safety valve
- Comprehensive alarms
- Low power consumption
- Low Maintenance





NMT-156
ACENTRIC OT TABLE

Your Search for truly OT Products Ends Here.

- FEATURES**
- Eccentric Pillars
 - 304 Stainless Steel Fitting
 - Detachable Head, Leg & Pelvic Section
 - Four Sectional Radio-translucent Top
 - Stand By Battery Backup With Re-charging Circuit.



Height Up	Head Up	Tender Burg	Tilt
Height Down	Head Down	Reverse Tenderburg	

NMT-157
CENTRIC OT TABLE

Your Search for truly OT Products Ends Here.

- Electric Gynecology And Obstretic OT Table





NMT-158
ELECTRO SURGICAL UNIT (Cautery Machine)

Your Search for truly OT Products Ends Here.

Features

- Microprocessor Based
- Smooth And Pure Cut
- Good Hemostasis Effect
- Minimum Cutting Effect
- Isolated Bipolar Output With No Sparking Characteristic
- Assured High Precision And Reliability In Surgical Procedures



- 400 Watt Digital Cautery Machine



- 300 watt Digital Cautery Machine
- 250 watt Digital Cautery Machine



- 300 Watt Analog Cautery Machine
- 250 Watt Analog Cautery Machine

NMT-159
SUCTION MACHINE

Your Search for truly OT Products Ends Here.

- Double Bottle Suction Machine (Indian)
Top SS (1/4HP)
- Double Bottle Suction Machine (Indian)
Top SS (1/2HP)
(Crompton Motor, Copper Motor)



- Double Bottle Suction Machine Imported



- Single Bottle Suction Machine Imported





CHAMBER DIMENSION	Automatic
20"x 48"(4 DRUM,248 Lit.)	Fully / Semi
20"x 36"(3 DRUM,180 Lit.)	Fully / Semi
16" x 24"(2 DRUM,81 Lit.)	Fully / Semi
18"x18"x36"(176 Liter)	Fully / Semi
24"x 24"x48"(450 Liter)	Fully

NMT-160
HORIZONTAL AUTO CLAVE

Your Search for truly OT Products Ends Here.
Manual Mode Also Available

- Chamber And All Contact Parts Made From SS 304 Grade Material
- PLC Based Automation
- Preprogrammed Air Removal Pocket Before Sterilization & Built Drying Facility After Sterilization To Get Your Material Dry.
- Central Locking System Provided For Safety During Machine In Operation Mode
- In Built Boiler Provided With Every Model.

Fully Automatic Features

- The Touch-screen For Visual Process Control And Connections For Measuring Signals
- Printers Automatic Printing Facility For Control In Case Of Heat Labile Materials
- Contamination Seal To Prevent Cross-flow Of Air Between The Clean And Unclean Area.

Optional

- Double Door With Sliding • Hinged Door
- Full Body Cover

NMT-161
VERTICAL AUTO CLAVE

Your Search for truly OT Products Ends Here.



- Single Drum Autoclave
Size : (12*11)
 - Double Drum Autoclave
Size : (12*22)
 - Triple Drum Autoclave
Size : (12*33)
- (Electric With Gas)**



Double Drum High pressure (Vertical)75 liter

- Triple Walled Unit Structure
- 304 Grade Stainless Steel
- 3KW Heter
- Automatic Low Water Cut Off
- Automatic Exhaust Value
- Working Pressure: 1.2 Kg / Cm2 & 2.1Kg / Cm2
- Temperature: 121°C & 134°C
- Double Drum High pressure
- Size : 16*24





NMT-162

ANESTHESIA TROLLEY

Your Search for truly OT Products Ends Here.

Anesthesia Trolley (1N2, 1O2)

Anesthesia Trolley (1N2, 2O2)

Anesthesia Trolley (2N2, 2O2)


- Anesthesia Trolley Accessory

Duo Top Circle Absorber

Single Top Circle Absorber

Close circuit port

Vaporizer Port



NMT-163

ANESTHESIA WORKSTATION

Your Search for truly OT Products Ends Here.

VCV, SIMV, ACMV, Manual, Standby,

Screen Size 15.5" Color LCD Display

VCV, PCV, SIMV, ACMV, Manual, Standby,

Screen Size 8.4" Color LCD Display

VCV, PCV, SIMV, ACMV, PSV, Manual,

Standby, Screen Size 12.1"

Color LCD Display

Halothane / Isoflurane /

Sevoflurane Vaporizer





NMT-164
MODULAR OT

Your Search for truly OT Products Ends Here.

OT Hight _____ OT Length _____ OT Width : _____

Distance Between AHU & Laminar : _____

NABH : YES / NO
Modular : _____
Semi Modular : _____



- Air Handling Unit & Condensing Unit

- 500CFM To 40000 CFM
Depending on room area
- Blower & Motor With Super Silent Centrifugal
- Blower Statically And Dynamically Balanced



- Ducting And Insulation Work
0.5 mm Duct In GI
(Make Jindal, Tata, Essar)
Insulation Of Air Duct With
9 mm Thick
Armorfoam Insulation.



- The Wall And Flooring
- PPGI/SS
 - Vinyl Flooring
 - Epoxy Flooring (Squire feet)
 - Corner Covering (Runing feet)

NMT-165
OT PENDENT

Your Search for truly OT Products Ends Here.

- Fix Arm OT Pendant - Single Arm OT Pendant - Laparoscopy Pendant
- Complete Alluminium



Feature • Tray-2 Nos • Switch-6 Nos • Oxygen-2 Nos • Vaccum-2 Nos • Air-4 Bar-1nos • Air-7 Bar-1 Nos





TM-100



TM-200

NMT-166

FETAL DOPPLER (TM-100 - TM-200)

Your Search for truly GYNEC Products Ends Here.

LCD Display Size 45mm.x25mm Ultrasound probe solid durable, anti-shock, anti-falling down Ultrasound frequency:2/2.5/3 MHz+3% Measuring Range - 50-240 BPM | Humidity:0-85% Ultrasound intensity:<2.5Mw/cm' Temperature:5c-40c, Rechargeable Ni-H battery (9v) Size - 135mm x 95mm x 35mm Weight -500gm.

NMT-167

WIRELESS FETAL TOUCH SCREEN MONITOR (ANDROID MODEL) (OX - 500 10.2 INCH)

- (1) High Resolution 10.2" TFT Screen, foldable touch screen display
- (2) Built-in 150mm thermal Printer
- (3) Super Sensitive wired or wireless probe, accurate and reliable results
- (4) Logn time data storage, can be played back and printed
- (5) Rechargeable battery support long time work
- (6) DSP technology, real-time identifying and measurement of FHR
- (7) Automatic Fetal Movement Dectection and Accounting
- (8) CTG Scoring system : Krebs Score/Fischer Score
- (9) Wire/Wireless connection with Central Monitoring System



NMT-168

FETAL MONITOR (BISTOS-BT-350)

Your Search for truly NICU Products Ends Here.

- (1) 7" wide TFT color LCD Display
- (2) Up to 150 patients data saving
- (3) Fetal heart sound play and record in PC
- (4) Rotation Screen
- (5) Waterproof Probe
- (6) GE/Philips Paper Compatible
- (7) Interpretation Function
- (8) Multiple language support
- (9) Built-in quick guide (LCD Type) Waterproof Probe
- (10) Twins Probe Available





ICU BED (MOTORIZED)



ICU BED (MANUALLY)



FULL FOWLER HI-LOW BED



FULL FOWLER BED



SEMI FOWLER HI-LOW BED



SEMI FOWLER BED



BACK-REST BED



GENERAL BED



GENERAL BED WITH M.S. RAILING



GENERAL BED WITH LAMINATED PANEL



CRIB CARRIER BED



DELUXE LOCKER





MONITOR TROLLEY



INTRUMENT TROLLEY (S.S)



TROLLEY WITH BOWL (S.S)



ECG TROLLEY (S.S)



TWO SECTION LABOR TABLE (S.S)



LABOR TABLE



EXAMINATION CUM GYNAEC TABLE



EXAMINATION TABLE



CARDIAC TABLE WITH HANDLE SYSTEM



FOLDING STRETCHER TROLLEY



SOILED LINEN TROLLEY (S.S)



VISITOR SEAT CUM BED WITH MATTRESS





Full Face Mask



Nasal Mask



Nasal Pillow



Tube



Wireless Spo2



Heating Tube



TPU Cover Nibp Cuff



Nylon Cover NIBP CUFF



Sing Tube Double Tube





OXYGEN SENSORS



3PIN MOLEX



DRAGER



STEREO JACK



JUMBO



MAQUET



SLE

FLOW SENSORS



HAMILTON



DRAGER



GE

VENTILATOR CIRCUITS



REUSABLE CIRCUIT



DISPOSABLE CIRCUITS



BAIN CIRCUITS





PATIENT MONITOR ACCESSORIES



SP02 SENSORS



SP02 CONNECTORS



SP02 EXTENSION CABLE



NIBP CONNECTORS



NIBP HOSEPIPE



ECG CABLE



ECG LEADSET



ECG CABLE CLIP



ECG ELECTRODES

ECG CABLES



REUSABLE CIRCUIT





TEMPERATURE SENSORS



ESU ACCESSORIES



REUSABLE ELECTROSURGICAL

MONOPOLAR HANDLE WITH CABLE

ELECTROSURGICAL PENCIL



BIPOLAR PAD CABLE

BIPOLAR CABLE VALLEYLAB

BIPOLAR CABLE VALLEYLAB EUROPEAN



DISPOSABLE BIPOLAR CABLE

BIPOLAR FORCEPS

BIPOLAR FORCEPS FOR





All Type Hospital Equipment Rental Services - Medical equipment rental providers are not held to the same industry standards as OEMs. Since quality costs money, you can imagine the corners that get cut. As a result, patient care suffers, and your staff grows more frustrated every day. Quality and reliability problems also lead to “bad rental” habits that increase your total costs.





All Type Hospital Machines Repairing Services - We are a renowned Medical Equipment Repair Service Provider sited in Gujarat, India. Our offered Medical Equipment Repairing Services comprises support in Monitor, Syringe Pump, and Infusion Pump, Bipap, Siemens Ventilators, Dragger Ventilators, Defibrillator, OT Table, OT Light, Laboratory equipment. The company is equipped with the latest tool and techniques, and provide its clients with a component level repairing, which is very cheaper than PCB level repairing. Clients can approach us for obtaining our Medical Equipment Repairing Services, at fair prices.



Contact Us



Office Address

Tirupati Complex, FF-108, Sardar Chowk, opp.
Ankur hospital, Krishnanagar, Ahmedabad,
Gujarat 380038



Phone Number

+91-8460675102



Email Id

nationalmt17@gmail.com

